



Bonding and banding products



Global Experts in Orthodontics

03

GC Orthodontics
GLOBAL EXPERTS IN ORTHODONTICS

Bonding & banding Summary

GC Ortho Connect	152
GC Ortho Etching Gel	152
GC Kommonbase	153



Produits de Collage et de scellement

Ces produits de collage sont innovants et apporteront un plus dans votre pratique quotidienne. Elaborés par GC au Japon, ils ont fait l'objet de tests approfondis afin de garantir une bonne tenue de vos attaches à coller en métal, céramique ou plastique.



Light cure adhesive



GC Ortho Connect

One step system

No need to apply primer or adhesive on the tooth, as it is already incorporated into the paste. A simple etching of the enamel followed by drying is enough to prepare the bracket or buccal tube for bonding.

	P/N
2 syringes of 2.9 g, 4 dispensing tips, 2 light protective covers	90-1100-0000



GC Ortho Etching Gel

With an ideal fluidity, this etching in gel (Phosphoric acid) will have a perfect consistency during the necessary time to prepare the tooth enamel.

	P/N
2 syringes of 4,8 g and 40 dispensing tips	90-1400-0000



Orthodontics lingual base extensions

GC Kommonbase



Improvements in customized base design and bonding materials have led to the development of a precise direct-bonding system called GC KommonBase. A large bonding base provides more precise bracket fit and positioning and maximizes bond strength. No transfer trays are needed because of the base's self-positioning shape. In addition, by using a reference jig, single brackets can be rebonded in case of failure or the need for repositioning.

2 syringes of 2,7 g, 4 dispensing tips and 1 light protective cover	P/N
KommonBase HV Clear (high viscosity)	90-1200-0000
KommonBase LV Pink (low viscosity)	90-1300-0000

